

# ADVENTURER PRECISION – AX1502 – 30127966

## SPECIFICATION

Maximum Capacity	1,520 g
Readability	0.01 g
Pan Size	7.7 in x 6.9 in (195 mm x 175 mm)
Internal Calibration	AutoCal™ - Automatic
Draftshield	—
Auto Door	Not Applicable
Auxiliary Display Model	Available as an Accessory
Battery Life	Not Applicable

---

Communication	USB; RS232 (Included); USB Host
Dimensions	13.93 in x 3.94 in x 9.05 in (354 mm x 100 mm x 230 mm) (LxHxW)
Display	High resolution full colour VGA touch display
In-Use Cover	Included
Legal for Trade	Not Applicable
Linearity ±	0.02 g
Minimum Weight (USP, 0.1%, typical)	20 g
Net Weight	10.2 lb (4.6 kg)
Pan Construction	Stainless Steel

---

---

Power	AC Adapter (Included)
Repeatability, typical	0.01 g
Stabilization Time	1.5 s
Tare Range	To capacity by subtraction
Units of Measurement	Singapore Tael; Ounce Troy; Pennyweight; Grain; Tical; Custom; Milligram; Newton; Momme; Taiwan Tael; Baht; Gram; Hong Kong Tael; Pound; Tola; Mesghal; Carat; Ounce
Working Environment	50°F – 85°F, 80%RH, non-condensing (10°C – 30°C, 80%RH, non-condensing)

---