

GUARDIAN 2000 HOTPLATE & STIRRER E-G21HS07C - 30680270

SPECIFICATIONS

Function	Heating and Stirring
Capacity	15 L
Plate Construction	Ceramic
Plate Size	7 in x 7 in (178 mm x 178 mm)
Temp. Range	70°C - 500°C
Speed Range	200 rpm – 1,600 rpm
Control	Analog
Dimensions	12.5 in x 4.4 in x 8.8 in (318 mm x 112 mm x 223 mm) (LxHxW)

Power Consumption	1,035 W
Working Environment	5°C – 40°C, 80% RH, Non-condensing (41°F – 104°F, 80% RH, Non-condensing)
Speed Range	200 to 1600 rpm
Net Weight	8.6 lb (3.9 kg)
Safety Certification	TUV
Power	220V-240V~ 7A 50Hz

