

# VORTEX MINI MIXERS - VNMNAL – 30392117

## SPECIFICATIONS

Control	Analog
Duty Rating	Intermittent
Orbit	0.2 in (4.9 mm)
Speed Range	300 rpm – 2800 rpm
Capacity	0.5 kg
Timer	Not Applicable
Dimensions	8.3 in x 6.5 in x 4.8 in (21 cm x 16.5 cm x 12.3 cm) (LxHxW)
Net Weight	9.7 lb (4.4 kg)

---

Power	230V, 0.6A, 50/60Hz
Power Consumption	135 W
Safety Certification	CE; TUV
Working Environment	39°F – 104°F, 20% – 85%RH, non-condensing (4°C – 40°C, 20% – 85%RH, non-condensing)

---